

# Human iPSC-Derived Cardiac Organoid differentiation Kit

## Human iPSC-Derived Cardiac Organoid Differentiation Kit

Cat. No. : RIPO-HWM002K

### Product Description

Cardiac organoids are three-dimensional in vitro models with a cellular composition and structural organization that is representative of the human heart. Human iPSC-Derived Cardiac Organoid maintenance Kit allows the maturation and the long-term maintenance of cardiac organoids.

### Product Specification

The basic medium of this differentiation kit is a serum-free, well-defined medium with minimal batch variation to which differentiation factors are added. This medium does not contain antibiotics, the addition of which may affect cardiac organoid differentiation.

### Product Information

Name	Component #	Size	Storage	Shelf Life
Basal Medium M-M	RIPO-HWM004-C01	225ml	4 °C	Stable for 1 years from date of manufacture (MFG) on label
Supplement M-M	RIPO-HWM004-1-C01	25ml	-20 °C	Stable for 1 years from date of manufacture (MFG) on label


### Materials Required but Not Included

- Ultra-Low Attachment 96 Well Plate (Corning, #7007)
- Ultra-Low Attachment 6 Well Plate (Corning, #3471)
- Orbital shaker (any brand, 2 cm shaking diameter)

### Equipment Required

- Incubator (37°C, 5% CO<sub>2</sub>)
- Low-speed centrifuge with a swinging bucket rotor with an adaptor for plate holders
- Orbital shaker
- Biosafety cabinet

### Protocol Diagram



*Figure 1. Cardiac Organoid Differentiation Process*

The color differs each component of differentiation kit. The dashed line represents the time for medium changes. Morphology of cardiac organoid at each stage of differentiation could be observed.

## Preparation of Media

Use sterile technique when performing the following manipulation

Medium	Component	Volume	IN-USE STORAGE/STABILITY
Medium M-M (500ml)	Basal Medium M-M	225ml	Mix completely the Basal Medium M-M and Supplement M-M to get Medium M-M. Store at 2 - 8°C for up to 2 weeks or aliquot as desired.
	Supplement M-M	25ml	

*Note: Please do not heat the complete medium (mixture of basal medium and supplement). Use it directly as cold as 2-8 °C.*

## Directions for Use

Please read the entire protocol before proceeding.

Use sterile technique when performing the following protocols.

**Note: This kit only serves for the Maturation and Maintenance of cardiac organoids (refer to step 6 and 7 of the Cardiac Organoid Differentiation part). For the differentiation of cardiac organoid, please use Human iPSC-Derived Cardiac Organoid Differentiation Kit (RIPO-HWM002K).**

**Note: Before cardiac organoid culturing, please make sure that the culture system you use is mTeSR-based. If your culture system is not mTeSR, please make sure that you have transferred your cells to the mTeSR system for at least 4 passages.**

### EB Formation

1. Aspirate medium from hPSC culture and wash the well with 3 mL of pre-warmed D-PBS (Without Ca<sup>++</sup> and Mg<sup>++</sup>) 3 times.
2. Aspirate PBS and add 2 mL of Gentle Cell Dissociation Reagent.
3. Incubate about 7-10 minutes for digestion of iPSCs to single cells.

*Note: Incubation time may vary when using different cell lines or different cell dissociation reagents.*

4. Using pipettes to pipet cells for obtaining single cells and centrifuge at 300g, 4 °C for 5 minutes
5. Remove the supernatant and add 1 ml mTeSR Plus to resuspend cells.
6. Count cells using Trypan Blue and a hemocytometer.
7. Transfer appropriate number of cells into Medium F to acquire final concentration of 7,500 cells/100 µL.
8. Add 100 µL of cell suspension into each well of a 96-well round-bottom ultra-low attachment plate.
9. Centrifuge the ultra-low attachment plate seeded with cells at 300g, 4 °C for 5 minutes.
10. Incubate the plate at 37°C, 5% CO<sub>2</sub> for 24h to formation of embryoid bodies.
11. Observe plate under microscope. The size of formed EBs should be ranging from 300 to 500 µm and with a smooth round edge (see protocol diagram). Count this day as day 0.

## EB Validation


### Morphology

- Round in shape with smooth edges and with limited dead cells surrounding.
- Diameter should be between 400-500  $\mu\text{m}$ .



## Cardiac Organoid Differentiation

1. At day 0, remove all the medium in the wells. Add 200 $\mu\text{l}$  of Medium A at each well and incubate at 37°C, 5% CO<sub>2</sub> for 48h.
2. After 48h, remove the 200 $\mu\text{l}$  medium in each, add 200 $\mu\text{l}$  of medium B in each well and incubate at 37°C, 5% CO<sub>2</sub> for 24h.
3. Repeat step (2) 3 times (add 4 times of medium B in total)
4. After the last incubation with medium B, remove the 200 $\mu\text{l}$  of medium B, add 200 $\mu\text{l}$  of medium C and incubate at 37°C, 5% CO<sub>2</sub> for 24h.
5. Repeat step (4) 4 times (add 5 times of medium C in total)
6. After the last incubation with medium C, transfer all cardiac organoids into ultra-low attachment 6 well plate (the maximum number is 24 organoids per well) and add 5 ml **medium M-M** per well. Then put the plate on an orbital shaker (as shown figures), which was placed inside the incubator, with the speed of 100 rpm.




7. Change the **medium M-M** fully every other day with the volume of 5 ml.

## Differentiation Validation

### Directly after stage A:

- The edge of EBs would change to irregular or smooth (Both could be accepted).
- Size is ranging from 400 to 500  $\mu\text{m}$ .



Stage B:

- After first 24 h incubation with medium B, the edge would change to smooth in shape.
- After 72 h incubation with medium B, the EBs would appear whole vacuum in the middle (If the shape changes to irregular and no vacuum appears, the differentiation is failed).
- The size after B stage would be larger than 800  $\mu\text{m}$ .




Beating:


- Regular beating would appear at 6- 12 days from EB formation.
- If no beating is observed even day 12 from differentiation, the differentiation could be considered as failed.

Marker expression (for matured cardiac organoids)


- Presence of cardiomyocytes: cTNT marker expression, which occupies almost 80 % of total cells, is acceptable (Day 15).



- Presence of endothelia cell: CD31 (Day 25)




- Presence of ventricle and atrium chambers: MLC2V and MLC2A respectively (Day 10)



## Related Products


Product	Cat. No.
Cardiac Organoid maintenance medium	RIPO-HWM002

## Validation Data of Cardiac Organoids




*Figure 2. Immunostaining of Cardiac Organoids*

(A) 20x view of cardiac organoids with cardiomyocytes were visualized using the cTNT marker with nuclei visualized by DAPI. (B) Presence of cardiomyocytes throughout the entire organoid can be observed.



*Figure 3. Transient Voltage Potential across Cardiac Organoids*

Silicon probes were placed across a cardiac organoid to measure the impulse generation and propagation causing the contraction activity. Each contraction lasts around 449 ms and can be continuously observed.



*Figure 4. AAV Capsid Screening by Cardiac Organoids*

Cardiac organoids grown for 11 days were infected with (A) AAV5-WT, (B) IVB-1, (C) IVB-2. Each of them was individually placed into a well plate. Transgene delivery efficacy was visualized by fluorescent intensity from GFP transgene expression.