Catalog # GMP-L18H16



#### Features

- Designed under ISO 9001:2015 and ISO 13485:2016
- Manufactured and QC tested under a GMP compliance factory
- Animal-Free materials
- Beta-lactam materials free
- Batch-to-batch consistency
- Stringent quality control tests
- No animal derived peptone and lactose used in production process

#### Source

GMP Human IL-18 Protein(GMP-L18H16) is expressed from E. coli cells. It contains AA Tyr 37 - Asp 193 (Accession # <u>Q14116</u>). Predicted N-terminus: Tyr 37

## **Molecular Characterization**

IL-18(Tyr 37 - Asp 193) Q14116

This protein carries no "tag".

The protein has a calculated MW of 18.2 kDa. The protein migrates as 17 kDa±3 kDa when calibrated against <u>Star Ribbon Pre-stained Protein Marker</u> under reducing (R) condition (SDS-PAGE).

# Endotoxin

Less than 10 EU/mg by the LAL method.

## **Host Cell Protein**

<0.5 ng/µg of protein tested by ELISA.

## **Host Cell DNA**

<0.02 ng/µg of protein tested by qPCR.

## **SDS-PAGE**



## Sterility

The sterility testing was performed by membrane filtration method described in CP<1101>, USP<71> and Eur. Ph. 2.6.1.

## Mycoplasma

Negative.

Purity

>95% as determined by SDS-PAGE.

## Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with protectants.

Contact us for customized product form or formulation.

#### Shipping

*This product is supplied and shipped with blue ice, please inquire the shipping cost.* 

## Storage

Upon receipt, store it immediately at -20°C or lower for long term storage.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 5 years in lyophilized state;
- -70°C for 12 months under sterile conditions after reconstitution.








#### Catalog # GMP-L18H16

GMP Human IL-18 Protein on SDS-PAGE under reducing (R) and nonreducing (NR) conditions. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With <u>Star Ribbon Pre-stained Protein</u> <u>Marker</u>).


# **Bioactivity-Bioactivity CELL BASE**

GMP Human IL-18 Protein stimulates secretion of IFN-γ by KG-1 cells.




GMP Human IL-18 Protein (Cat. No. GMP-L18H16) stimulates secretion of IFN- $\gamma$  by KG-1 cells. The specific activity of GMP Human IL-18 Protein is > 3.00x10^6 IU/mg, which is calibrated against WHO Reference Reagent Interleukin-18 (Human rDNA derived) (NIBSC code: 03/200) (QC tested).

# **Application Data**



Human PBMCs were cultured with GMP Human IL-2 Protein (ACROBiosystems, Cat. No. GMP-L02H14), GMP Human IL-18 Protein (ACROBiosystems, Cat. No. GMP-L18H16) and GMP Human 4-1BB Ligand Protein (ACROBiosystems, Cat. No. GMP-41LH26), in CelThreaTM GMP T Cell Expansion Medium (ACROBiosystems, Cat. No. GMP-CM3101 & GMP-CM3101-1) for two weeks. The result shows that GMP Human IL-18 Protein can promote the expansion of these cells with a reasonable cell viability. Notably, the cells exhibit better expansion in GMP Human IL-18 Protein (ACROBiosystems).

GMP Human IL-18 Protein stimulates the secretion of IFN- $\gamma$  by KG-1



The activity of GMP Human IL-18 Protein (Cat. No. GMP-L18H16) was

higher than other competing products.




>>> www.acrobiosystems.com




Catalog # GMP-L18H16





Human PBMCs were cultured with GMP Human IL-2 Protein (ACROBiosystems, Cat. No. GMP-L02H14), GMP Human IL-18 Protein (ACROBiosystems, Cat. No. GMP-L18H16) and GMP Human 4-1BB Ligand Protein (ACROBiosystems, Cat. No. GMP-41LH26), in CelThreaTM GMP T Cell Expansion Medium (ACROBiosystems, Cat. No. GMP-CM3101 & GMP-CM3101-1) for two weeks. The result shows that GMP Human IL-18 Protein (ACROBiosystems) can be comparable with the Competitor M IL-18 protein.

# **Bioactivity-Stability**



The Cell based assay shows batch-to-batch consistency between Acro's GMP and PG IL-18.

# MANUFACTURING SPECIFICATIONS

ACROBiosystems GMP grade products are produced under a quality management system and in compliance with relevant guidelines: Ph. Eur General Chapter 5.2.12 Raw materials of biological origin for the production of cell-based and gene therapy medicinal products; USP<92>Growth Factors and Cytokines Used in Cell Therapy Manufacturing; USP<1043>Ancillary Materials for Cell, Gene, and Tissue-Engineered Products; ISO/TS 20399-1:2018, Biotechnology - Ancillary Materials Present During the Production of Cellular Therapeutic Products.

ACROBiosystems Quality Management System Contents:

Designed under ISO 9001:2015 and ISO 13485:2016, Manufactured and QC tested under a GMP compliance factory.

Animal-Free materials

Materials purchased from the approved suppliers by QA

ISO 5 clean rooms and automatic filling equipment



>>> www.acrobiosystems.com

10/12/2024

Catalog # GMP-L18H16

#### Qualified personnel

Quality-related documents review and approve by QA Fully batch production and control records Equipment maintenance and calibration Validation of analytical procedures Stability studies conducted Comprehensive regulatory support files

Request For Regulatory Support Files (RSF)

ACROBiosystems provide rigorous quality control tests (fully validated equipment, processes and test methods) on our GMP grade products to ensure that they meet stringent standards in terms of purity, safety, activity and inter-batch stability, and each bulk QC lot mainly contains the following specific information:

# Protein content Endotoxin level Residual Host Cell DNA content Residual Host Cell Protein content Biological activity analysis Microbial testing

Mycoplasma testing

**SDS-PAGE** 

In vitro virus assay

Residual moisture

Batch-to-batch consistency

## Background

Interleukin-18 (IL-18) is a potent proinflammatory cytokine that induces interferon-gamma (IFN-gamma) production from Th1 cells, NK cells and activated macrophages, particularly in the presence of IL-12. IL-18 also functions in developmental regulation of T-lymphocyte helper type I cells and in Fas-mediated cytotoxicity. Suppression of IL-18 activity is being investigated for treatment of chronic inflammatory diseases such as Crohn's disease and rheumatoid arthritis. It acts by inducing heterodimerization of the two subunits of its receptor, IL-18RAlpha and IL-18RBeta shows structural similarity to IL-1.

## **Clinical and Translational Updates**





>>> www.acrobiosystems.com

