

# Human ADAM9 Protein, His Tag (active enzyme)

Catalog # AD9-H52H7



BIOSYSTEMS  
**Acro**  
Surprise Inside!

## Synonym

Disintegrin and metalloproteinase domain-containing protein 9 (EC:3.4.24.-)  
Cellular disintegrin-related protein, Meltrin-gamma, Metalloprotease, disintegrin, cysteine-rich protein 9, Myeloma cell metalloproteinase, ADAM9, KIAA0021, MCMP, MDC9, MLTNG

## Source

Human ADAM9, His Tag(AD9-H52H7) is expressed from human 293 cells (HEK293). It contains AA Ala 206 - Asp 697 (Accession # [Q13443-1](#)).  
Predicted N-terminus: Ala 206

## Molecular Characterization

ADAM9(Ala 206 - Asp 697)  
Q13443-1 Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 55.2 kDa. The protein migrates as 65-70 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

## Endotoxin

Less than 1.0 EU per µg by the LAL method.

## Purity

>95% as determined by SDS-PAGE.

## Formulation

Supplied as 0.2 µm filtered solution in 20 mM Tris, 500 mM NaCl, pH7.5 with glycerol as protectant.

Contact us for customized product form or formulation.

## Shipping

*This product is supplied and shipped with dry ice, please inquire the shipping cost.*


## Storage

*Please avoid repeated freeze-thaw cycles.*

This product is stable after storage at:

- The product MUST be stored at -70°C or lower upon receipt;
- -70°C for 3 months under sterile conditions.

## SDS-PAGE



Human ADAM9, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With [Star Ribbon Pre-stained Protein Marker](#)).

## Bioactivity-ELISA

Discounts, Gifts,  
and more!



➤ [www.acrobiosystems.com](http://www.acrobiosystems.com)



5/30/2024

# Human ADAM9 Protein, His Tag (active enzyme)

Catalog # AD9-H52H7



BIOSYSTEMS  
**Acro**



Immobilized Human ADAM9, His Tag (Cat. No. AD9-H52H7) at 1 µg/mL (100 µL/well) can bind Izeltabart with a linear range of 0.06-1 ng/mL (QC tested).

Immobilized Izeltabart at 1 µg/mL (100 µL/well) can bind Human ADAM9, His Tag (Cat. No. AD9-H52H7) with a linear range of 2-31 ng/mL (Routinely tested).

## Bioactivity

Measured by its ability to cleave a fluorogenic peptide substrate Mca-PLAQAV-Dpa-RSSSR-NH<sub>2</sub>. The specific activity is >10 pmol/min/µg (QC tested).

## Background

ADAM9 (A disintegrin and a metalloprotease 9) is a membrane-anchored protein that participates in a variety of physiological functions, primarily through the disintegrin domain for adhesion and the metalloprotease domain for ectodomain shedding of a wide variety of cell surface proteins. ADAM9 influences the developmental process, inflammation, and degenerative diseases. Recently, increasing evidence has shown that ADAM9 plays an important role in tumor biology. Overexpression of ADAM9 has been found in several cancer types and is correlated with tumor aggressiveness and poor prognosis. In addition, through either proteolytic or non-proteolytic pathways, ADAM9 promotes tumor progression, therapeutic resistance, and metastasis of cancers. Therefore, comprehensively understanding the mechanism of ADAM9 is crucial for the development of therapeutic anti-cancer strategies. In this review, we summarize the current understanding of ADAM9 in biological function, pathophysiological diseases, and various cancers. Recent advances in therapeutic strategies using ADAM9-related pathways are presented as well.

## Clinical and Translational Updates

Discounts, Gifts,  
and more!

